



ResMed

Actual Size: 5 1/2" (Width)
x 6 7/8" (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirFit™ F20

Nasal mask
Nasal mask for Her

© 2016 ResMed Ltd. All Rights Reserved.
Specifications may change without notice.
Product availability may vary across regions.
1031822 / 1 2016-09

CE
Rx Only 0123

Full face mask
Full face mask for Her

AirFit™ F20

SIZE
S M L

0 Inch
1
2
3
4
5
6
7

© 2016 ResMed Ltd. All Rights Reserved. Specifications may change without notice.
Product availability may vary across regions.

ResMed

Rx Only **CE** 0123